



- 材料是有机硅和芦荟相结合，因此内部舒适度优于纯硅胶材质。
- 适用于皮肤敏感的患者使用，例如植皮、有骨突或易过敏的皮肤。
- 外部织物提供支撑和保护，同时保持柔韧性和弹性。
- 可以在 30° C 水温下手洗或机洗。
- 分为带锁孔和不带锁孔两种，其中带锁孔的硅胶套下端织物有 10C M 的网状交叉织物增强。



残肢末端向上 4cm
测量围长确定尺寸





带锁孔

型号	厚度	末端厚度	围长
1S711-3-160	3 mm	17,5 to 18,5 mm	16
1S711-3-180			18
1S711-3-200			20
1S711-3-220			22
1S711-3-235			23.5
1S711-3-250			25
1S711-3-265			26.5
1S711-3-280			28
1S711-3-300			30
1S711-3-320			32
1S711-3-340			34
1S711-3-360			36
1S711-3-400			40
1S711-3-450			45
1S711-6-160			6 mm
1S711-6-180	18		
1S711-6-200	20		
1S711-6-220	22		
1S711-6-235	23.5		
1S711-6-250	25		
1S711-6-265	26.5		
1S711-6-280	28		
1S711-6-300	30		
1S711-6-320	32		
1S711-6-340	34		
1S711-6-360	36		
1S711-6-400	40		
1S711-6-450	45		

不带锁孔

型号	厚度	末端厚度	围长		
1S721-3-160	3 mm	14,5 to 16 mm	16		
1S721-3-180			18		
1S721-3-200			20		
1S721-3-220			22		
1S721-3-235			23.5		
1S721-3-250			25		
1S721-3-265			26.5		
1S721-3-280			28		
1S721-3-300			30		
1S721-3-320			32		
1S721-3-340			34		
1S721-3-360			36		
1S721-3-400			40		
1S721-6-160			6 mm	16 to 17,5 mm	16
1S721-6-180					18
1S721-6-200	20				
1S721-6-220	22				
1S721-6-235	23.5				
1S721-6-250	25				
1S721-6-265	26.5				
1S721-6-280	28				
1S721-6-300	30				
1S721-6-320	32				
1S721-6-340	34				
1S721-6-360	36				

- 材料是有机硅。
- 适用于皮肤敏感的患者使用，例如植皮、有骨突或易过敏的皮肤。
- 外部织物提供支撑和保护，同时保持柔韧性和弹性。
- 可以在 30° C 水温下手洗或机洗。
- 分为带锁孔和不带锁孔两种，其中带锁孔的硅胶套下端织物有 10CM 的网状交叉织物增强。



残肢末端向上 4cm
测量围长确定尺寸



带锁孔

型号	厚度	末端厚度	围长
1S712-3-160	3 mm	17,5 to 18,5 mm	16
1S712-3-180			18
1S712-3-200			20
1S712-3-220			22
1S712-3-235			23.5
1S712-3-250			25
1S712-3-265			26.5
1S712-3-280			28
1S712-3-300			30
1S712-3-320			32
1S712-3-340			34
1S712-3-360			36
1S712-3-400			40
1S712-3-450			45
1S712-6-200			6 mm
1S712-6-220	22		
1S712-6-235	23.5		
1S712-6-250	25		
1S712-6-265	26.5		
1S712-6-280	28		
1S712-6-300	30		
1S712-6-320	32		

不带锁孔

型号	厚度	末端厚度	围长
1S722-3-160	3 mm	14,5 to 16 mm	16
1S722-3-180			18
1S722-3-200			20
1S722-3-220			22
1S722-3-235			23.5
1S722-3-250			25
1S722-3-265			26.5
1S722-3-280			28
1S722-3-300			30
1S722-3-320			32
1S722-3-340			34
1S722-3-360			36
1S722-6-160	6 mm	16 to 17,5 mm	16
1S722-6-180			18
1S722-6-200			20
1S722-6-220			22
1S722-6-235			23.5
1S722-6-250			25
1S722-6-265			26.5
1S722-6-280			28
1S722-6-300			30
1S722-6-320			32
1S722-6-340	34		
1S722-6-360	36		